5255462011 OCF DE

Last updated 2020-11-12 Rev.B

Resin-based Dental Restorative Material

OMNICHROMA *Flow*

CE ¹⁶³⁹ REF OC02

SYRINGE

Instructions for Use

Manufacturer

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

EC Rep.

Tokuyama Dental Italy S.r.l.

Via Chizzalunga, 1

36066 Sandrigo, Vicenza, Italy

Tel: +39 0444 659650

Fax: +39 0444 750345

Importer

Tokuyama Dental Deutschland GmbH

Am Landwehrbach 5

48341 Altenberge Germany

Tel: +49 2505 938513

Fax: +49 2505 938515

ENGLISH	2	ČESKÝ	7	DANSK	13	DEUTSCH	20
EESTI KEEL	27	FRANÇAIS	33	ITALIANO	39	LATVISKI	46
LIETUVIŲ K.	52	MAGYAR	57	NEDERLANDS	64	NORSK	70
POLSKI	75	SLOVENČINA	82	SUOMI	88	SVENSKA	94

ENGLISH

Carefully read all information, precautions, notes and materials before using.

PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) OMNICHROMA FLOW is a low viscosity, light-cured, radiopaque composite resin for use in anterior and posterior restorations and is indicated for all carious classes. OMNICHROMA FLOW is a single shade material. Direct placement of OMNICHROMA FLOW into a prepared cavity allows for easy handling.
- 2) OMNICHROMA FLOW contains 71% by weight (57% by volume) of spherical silica-zirconia filler (mean particle size: 0.3µm, particle size range: 0.2 to 0.4µm) and composite filler. OMNICHROMA FLOW contains 1,6-bis(methacryl-ethylloxycarbonylamino)trimethyl hexane (UDMA), 1,9-nonamethylene glycol dimethacrylate, Mequinol, Dibutyl hydroxyl toluene and UV absorber.
- 3) OMNICHROMA FLOW is available in SYRINGE.

INDICATIONS

- Direct anterior and posterior restorations
- Cavity base or liner
- Repair of porcelain/composite

CONTRAINDICATIONS

OMNICHROMA FLOW contains methacrylic monomers and UV absorber. DO NOT use OMNICHROMA FLOW for patients allergic to or hypersensitive to methacrylic and related monomers, UV absorber or any of the other ingredients.

PRECAUTIONS

- 1) DO NOT use OMNICHROMA FLOW for any purpose other than those listed in these instructions. Use OMNICHROMA FLOW only as directed herein.
- 2) OMNICHROMA FLOW is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- 3) DO NOT use OMNICHROMA FLOW if the safety seals are broken or appear to have been tampered with.
- 4) If OMNICHROMA FLOW causes an allergic reaction or hypersensitivity, discontinue its use immediately.

- 5) Use examination gloves (plastic, vinyl or latex) at all times when handling OMNICHROMA FLOW to avoid the possibility of allergic reaction from methacrylic monomers or UV absorber.

Note: Certain substance/materials may penetrate through examination gloves. If OMNICHROMA FLOW comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.

- 6) Avoid contact of OMNICHROMA FLOW with eyes, mucosal membrane, skin and clothing.
 - If OMNICHROMA FLOW comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If OMNICHROMA FLOW comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If OMNICHROMA FLOW comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) OMNICHROMA FLOW should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 8) To avoid the unintentional ingestion of OMNICHROMA FLOW, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments and brushes with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude. After removing from the refrigerator, allow the product to sit 15 minutes or until it reaches room temperature.
- 12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.

INDICATIONS FOR FILLING AND CURING

OMNICHROMA FLOW is designed to be cured by either a halogen or LED curing-light with a wavelength of 400-500 nm. Be sure to light-cure OMNICHROMA FLOW extra-orally and check the time needed for complete hardening of OMNICHROMA FLOW with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light type	Intensity (mW/cm ²)	Curing time (seconds)	Increment depth (mm) ¹⁾
Halogen	400	20	2.7
		40	3.8
	600	20	3.0
		800	2.9
LED	400	20	2.7
		40	3.7
	600	20	2.9
		10	2.8
	800	20	3.4
		1200	10
2300	6	3.1	

1) Increment depth was determined on the basis of test results performed in accordance with "depth of cure" of ISO 4049.

■ STORAGE

- 1) Store OMNICHROMA FLOW at temperatures between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use OMNICHROMA FLOW after the indicated date of expiration on the syringe package.

■ DISPOSAL

To safely dispose of excess OMNICHROMA FLOW, extrude unused portion from SYRINGE and light-cure before disposal.

■ CLINICAL PROCEDURE

Determine if this product is suitable for application to the case before use. The desirable esthetic outcome may not be obtained depending on the case (e.g. a case where very little tooth structure remains, the underlying tooth is extremely discolored, or the color of the underlying tooth is to be changed to another color, etc.).

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation

A rubber dam is the preferred method of isolation.

3. Cavity Preparation

Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II), bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- Scalloped bevel could be desirable depending on the esthetics and retention.

- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

4. Pulp Protection

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing OMNICHROMA FLOW. OMNICHROMA FLOW is not indicated for direct pulp capping.

5. Bonding System


Apply a bonding system according to its manufacturer's instructions.

- If using self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

- In case of porcelain/composite repairs, be sure to condition the surface of porcelain/composite with a silane coupling agent BEFORE performing the bonding procedure.

6. Dispensing

- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip (metal tip) onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tips are exclusively designed for the SYRINGE of OMNICHROMA FLOW. Use the attached Dispensing Tips (metal tips).
- Dispensing Tips (TOKUYAMA Dispensing Tip) are also sold separately.
- Use uncontaminated gloves when handling the syringe and attaching the Dispensing Tip.

7. Filling and contouring

- Cavities can be directly filled from the Dispensing Tip, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated curing depth at a time (please refer to the aforementioned table).
- DO NOT apply excessive force to extrude the material from the syringe immediately after removal from the refrigerator.
 - DO NOT mix OMNICHROMA FLOW with other brands of resin composite to avoid incomplete cure, entrapment of air bubbles or losing its color adjustment.
 - After extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.
 - This product is provided non-sterile and is not intended to be sterilized by the user.
 - Use OMNICHROMA BLOCKER or OMNICHROMA BLOCKER FLOW under OMNICHROMA FLOW for the direct restoration of extensive Class III or IV cavity. OMNICHROMA FLOW can appear semi-translucent when OMNICHROMA FLOW is not adequately surrounded by a cavity wall. If applied in direct restorations of extensive Class III or IV cavity, the cured OMNICHROMA FLOW may appear dark.
 - In the case of repair of porcelain / composite, the application of an opaquer to the fractured

surfaces as an initial layer is needed to hide a metal portion if it exists. Then, place OMNICHROMA FLOW onto the fractured surfaces.

- In the case of masking tooth discoloration, apply OMNICHROMA BLOCKER, OMNICHROMA BLOCKER FLOW or a dental color masking material such as an opaquer on the discolored surface to mask the discoloration, followed by the placement of OMNICHROMA FLOW.

8. Curing

- Light-cure each increment for at least the indicated time (see aforementioned table).
- If other brands of composite resins are layered over the cured composite, follow their instructions.

9. Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.

Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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