Oral & Microsurgery – Grasping & Holding Instruments

Needle Holders, Hemostatic Forceps, Tweezers & Bracket Holding Tweezers

Ray-Bright Industries
Sialkot Pakistan.

CMC MEDICAL DEVICES & DRUGS, S.L
C/ Horacio, Lengo n18, C. P 29006
Málaga-Spain
Phone+34 951 214 054
Email: info@cmcmmedicaldevices.com

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en-English

Attention

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

Intended Use:

These instruments are intended to be used for precise handling of tissues, sutures, and small objects in oral and microsurgical procedures. They enable secure grasping, controlled manipulation, and atraumatic handling of delicate structures.

The Instrument Specific Intended Use along with their variants is given below:

Needle Holders:

Intended to be used for securely grasping and guiding suturing needles during dental, oral, and surgical procedures, ensuring accurate wound closure and minimal tissue trauma.

Variants:

Mathieu (Straight & Curved): For general suturing in dental, oral, and microsurgical procedures

Mathieu-Olsen: Includes integrated scissors for suturing and cutting without instrument exchange.

Lichtenberg: For precise control in fine suturing procedures.

Castroviejo / Micro: For microsurgical and fine suture placement in delicate tissues.

Crile-Ryde: For firm grip on small to medium-sized needles.

Fine-Swedish: For handling very fine needles in microsurgical applications.

Mayo-Hegar: For handling larger suture needles in heavy tissue procedures.

Crile-Wood: For precise control when suturing with delicate or medium-sized needles.

Hemostatic Forceps:

Intended to be used for clamping blood vessels and tissues for effective hemostasis.

Variants:

Halstead-Mosquito:For small vessel clamping.

Classic Pean: For large vessel or thick tissue occlusion.

Allison: For delicate tissue handling.

Kelly: For medium vessel clamping.

Kocher: For heavy tissue grasping.

Steiglitz: For small fragment removal.

Adson Hemostatic: For small tissue control.

Backhaus: For draping and tissue grasping.

Tohoku: For microsurgical tissue handling.

Gemini: For deep vessel clamping.

Hartman: For fine blood vessel clamping

Tweezers & Bracket Holding Tweezers:

Tweezers are used for handling soft tissues, sutures, and small objects.

Variants:

Adson: For soft tissue grasping.

Flagg: For cotton roll handling.

Perry: For use in periodontal procedures.

Miller: For suture placement.

College: For general dental applications.

Surgical: For holding sutures or foreign objects.

Meriam Tweezers: For handling small materials.

Anatomic: For surgical tissue grasping

Taylor: For fine material manipulation.

Semken: For delicate sutures.

Crane-Kaplan: Used in periodontal procedures.

Debakey: For vascular and microsurgical procedures.

Gerald: For microsurgical grasping.

Microsurgical: For ultra-fine materials handling.

Endo-Tweezers: For endodontic material handling.

Bracket Holding Tweezers: For orthodontic bracket placement

 Our Oral & Microsurgery – Grasping & Holding Instruments are intended to be used by qualified dental/healthcare professionals only in a professional dental/healthcare facility.

Patient population:

These instruments are intended for use on adult and child patients undergoing oral and microsurgical procedures. The appropriate instrument type and size must be selected based on the patient's oral anatomy, the surgical site, and the specific procedural requirements.

For Available Configurations/Models please visit the Oral & Microsurgery section in our Catalogue

Features

Needle Holders:

- Secure suturing needle grip.
- Ratchet locking mechanism for stability.
- Straight and curved designs for access.
- Stainless steel for durability and corrosion resistance.
- Fine jaws for precision suturing.

Hemostatic Forceps:

- Multi-position locking for controlled clamping.
- Serrated jaws prevent slippage.
- Straight and curved variants.
- High-grade stainless-steel construction.
- Fine tips for minimal tissue trauma.

Tweezers & Bracket Holding Tweezers:

- Ultra-fine tips for precision.
- Ergonomic anti-slip grooves.
- Stainless steel for longevity.
- Angled and curved designs for accessibility.
- Some models feature spring-loaded/self-closing mechanisms.

Operating instructions

 These instruments are non-sterile and must be cleaned and sterilized before first use.

 Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

 Ensure that ratchet and hinge mechanisms are functioning properly before use. Avoid excessive force to prevent damage to the locking mechanism

For Needle Holders:

- Select the appropriate needle holder for the suture size and procedure.
- Position the needle securely within the jaws.
- Use a controlled grip for precise suturing.
- Engage the ratchet mechanism for stability.

For Hemostatic Forceps:

- Choose the appropriate forceps for the vessel size and tissue type.
- Clamp the vessel using the locking mechanism to control bleeding.
- Ensure proper engagement of the ratchet to maintain pressure.
- Release gently to prevent tissue damage

For Tweezers & Bracket Holding Tweezers:

- Use tweezers to grasp and position materials with precision.
- For bracket holding tweezers, ensure a firm grip on the orthodontic bracket before placement.
- Handle with care to prevent accidental slipping.

For Periodontal Knife:

- Ensure proper blade selection for the intended surgical procedure.
- Use gentle, controlled strokes to make precise incisions while minimizing tissue damage.
- Maintain a steady grip and correct angulation to ensure optimal tissue contouring.

Adverse Event

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

Cleaning & Sterilization

Cleaning

 Clean the device with care after each application to avoid the drying of contaminants.

Initial Cleaning

1. Remove any visible debris from the instrument using moisten a cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Use a soft brush and enzymatic cleaner to remove contaminants thoroughly.

Manual Cleaning

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

Automated Cleaning

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.
4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERILIZATION PROCEDURE.

 For needle holder and forceps, ensure that the hinge and ratchet mechanism is fully cleaned to prevent residue buildup that could affect movement.

 Ultrasonic cleaning may be performed in accordance with hospital or facility protocols. Ensure proper rinsing and drying after the process to prevent residue buildup.

Sterilization

Note: Ensure that the instruments are completely dry before sterilization.

Steam Autoclaving

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. Ensure that the hinges and ratchet mechanisms are lubricated with an appropriate medical-grade lubricant after sterilization to maintain smooth operation.
5. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
6. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.
7. For scissors and needle holders, periodically check the hinge movement to ensure continued smooth function.

Important Notes

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.
- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturer are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.
- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.

Disposal

Our Grasping & Holding Instruments for Oral and Microsurgeries: Needle Holders, Hemostatic Forceps, Tweezers & Bracket Holding Tweezers contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before disposition. Please follow the local and national regulations or healthcare facilities' defined disposal and waste management policies.

Disclaimer

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

Explanation of utilized symbols

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)
	Manufacturing Date
	Manufacturer
	Authorized representative in the European Community
	Importer

	Keep in a cool, dry place
	Caution, read instructions

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