X-Post®

EN	Instructions For Use	2
FR	Instructions d'utilisation	8
IT	Istruzioni per l'uso	14
DE	Gebrauchsanweisung	20
ES	Instrucciones de uso	26
PT	Instruções de utilização	32
NL	Gebruiksaanwijzing	38
SV	Användarinstruktioner	44
DA	Brugsanvisning	50
NO	Brukerinstruksjoner	56
LT	Naudojimo instrukcija	62
LV	Lietošanas norādījumi	67
ET	Kasutusjuhised	72
CS	Návod k použití	77
SK	Návod na použitie	83
HU	Használati utasítás	89
RO	Instrucțiuni de utilizare	95
SR	Uputstva za upotrebu	101
BS	Uputstva za upotrebu	107
CNR	Uputstva za upotrebu	113
SQ	Udhëzimet për përdorim	119
EL	Οδηγίες χρήσης	125
TR	Kullanım Kılavuzu	131
PT-BR	Instruções de uso	137
KO	사용 지침	143

X-Post®

FOR DENTAL USE ONLY

FOR PROFESSIONAL USE ONLY

NON-STERILE – MANDATORY PROCESSING BEFORE FIRST USE

SINGLE USE – REPROCESSING PROHIBITED

INSTRUCTIONS FOR USE

POSTS: SINGLE USE NON-STERILE IMPLANTABLE ENDODONTIC DEVICES

REF. 60667332, 60667333, 6066734, 6066735 (X-Post®)

1 PRODUCT DESCRIPTION

- Product description:
X-Post® are single use and non-sterile radiopaque, translucent fiber post. They are made of glass fibers embedded in an epoxy resin matrix. Their shape is cylindric-conical and are light-conducting. These devices have a color code on their upper extremity according to their dimensions. Fiber Posts are available in a full range of sizes (see section [4]).
- Clinical benefits:
The clinical benefit for the patient being treated with X-Post® is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
The products are indicated for treatment of endodontic disease.
- Intended purpose:
Posts are devices intended to retain the core material in a tooth with extensive loss of coronal structure.
- Intended patient population:
Patients who need root canal treatment.
- Intended user and environment of use:
Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

Insufficient residual dentin: at least 2 mm of tooth structure is required around the preparation.

4 DELIVERY FORM

The X-Post® brand comprises the following devices:

Posts

SKU	Product	Available sizes					Packaging and content
		N°X: color	Anchor		Head		
			Diameter (mm)	Length (mm)	Diameter (mm)	Length (mm)	
6066733X	X-Post®	N°1: Yellow	0.8	6.23	1.35	13.77	Blister with 5 devices
		N°2: Red	0.8	7.48	1.47	12.52	
		N°3: Blue	1.0	7.49	1.67	12.51	
		N°4: Green	1.0	9.29	1.83	10.71	

5 COMPOSITION

X-Post® fiber posts are made of glass fibers (81% ± 3% weight) embedded in an epoxy resin matrix (19% ± 3% weight).

6 COMPATIBLE DEVICES

- X-Post® Fiber Posts are used in combination primer adhesive combined with dual cure resin cement and following devices: Largo® Peeso Reamer and EasyPost® Precision Drill.
- Penetration drilling prior to root canal insertion of Posts shall be performed using Largo® Peeso Reamer REF. A0009 Nr.1 or Nr.2 instruments. Please refer to the Largo® Peeso Reamer brand IFU for the Instructions For Use related to these instruments.
- Precision drilling prior to root canal insertion of Posts shall be performed using EasyPost® - Precision Drill REF. C0601 instruments. Please refer to the EasyPost® brand IFU for the Instructions For Use related to these instruments.

7 PERFORMANCE CHARACTERISTICS

- Expected life of Posts: 10 years.
- Posts radio-opacity level: 200% Al equivalent (ISO 4049: Dentistry - Polymer-based restorative materials).

8 SAFETY NOTES

8.1 Warnings

- Strictly follow the Instructions For Use to minimize the following risks to the X-Post® devices, the patient and/or the user:
 - Breakage of instrument.
 - Contamination and cross-contamination.
 - Swallowing of working part of the instrument.
 - Post debonding.
 - Toxic or allergic reactions caused by processing residues.
- The use of inappropriate rotation speed, pressure and contact duration would speed up the wear of drills for posts, increase the risk of breakage and heat creation.

- Please refer to the Largo® Peeso Reamer brand IFU for the Warnings related to these instruments.

8.2 Precautions

General

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Processing step is mandatory before use according to section [10] for all the products mentioned in this IFU.

Posts

- Check the Posts devices before use, do not use and discard the devices if damaged.
- Shortening of the Posts devices should be done outside the mouth using a diamond-tipped disc.
- Do not touch the Posts devices with your fingers after cleaning / sterilization.

EasyPost® - Precision Drill instruments

- Check the Drills for Posts instruments before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Check the Drills for Posts instruments and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended purpose with the required safety level, instruments should be discarded.
- Before using any Drills for Posts instrument, make sure it is well connected to the contra-angle head.
- Always use minimal apical pressure. Never force the Drills for Posts instruments down the canal.

Largo® Peeso Reamer instruments

- Please refer to the Largo® Peeso Reamer brand IFU for the Precautions related to these instruments.

8.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

8.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

Keep away X-Post® devices from light and heat.

8.5 MRI (Magnetic Resonance Imaging) safety information

Fiber Posts are MR Safe. A patient with X-Post® devices can be safely scanned without specific conditions.

9 STEP BY STEP INSTRUCTIONS

- 1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
- 2) Select the correct X-Post® size, according to anatomical situation of tooth using the radiograph.
 - ⇒ At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
 - ⇒ Select the post size so that 2/3 of its length will be inserted in the root canal.
- 3) Select the appropriate Largo® Peeso Reamer and EasyPost® - Precision Drill sizes corresponding to the selected X-Post size as indicated in the following table:

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- 4) Remove the gutta percha and perform penetration drilling with the selected Largo® Peeso Reamer size.
 - ⇒ Drill as deep as necessary for the post to be inserted in the root canal at 2/3 of its length.
 - ⇒ At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length must be reduced.
- 5) Perform precision drilling with the selected EasyPost® - Precision Drill size.
- 6) Rinse the canal with water or air to remove drilling residues.
- 7) Check the selected X-Post® positioning in mouth.
- 8) Shorten the selected X-Post® outside the mouth to its final length using a diamond-tipped disc.
- 9) Immerse or soak the selected X-Post® with isopropanol for 30 seconds. Let the post dry for 1 minute.
- 10) For dentin conditioning prior to post cementation, refer to adhesive or cement manufacturer's instructions for use.
- 11) Insert the sealing adhesive or cement in the root canal and on the post itself following adhesive or cement manufacturer's instructions for use¹.
- 12) Insert the X-Post® in the root canal to 2/3 of its length.
- 13) Wait for the cement to harden, then remove any excess.

¹. The adhesion between the Posts devices and the Prime&Bond Active with SCA and Core-x Flow has been tested.

10 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- X-Post® devices are single use implantable devices, re-use and processing more than 1 time is prohibited. Re-use can increase the risk of cross contamination, poor cementation or breakage.

- Please follow the General Processing Instructions For Endodontic Products (dentsplysirona.com/ifu) for the disinfection, cleaning and sterilization steps.

11 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

12 EXPIRATION DATE

None.

13 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

14 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for processing. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

15 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information will be available at both the European database on medical devices (EUDAMED) and the IFU website, where it will be linked to the Basic UDI-DI.

EUDAMED

<https://ec.europa.eu/tools/eudamed>

IFU

dentsplysirona.com/ifu

The Basic UDI-DI numbers are:

- Posts Fibers: ++J00310056E8

16 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

Authorised representatives		Manufacturer
 2797	EU REP DENTSPLY DeTrey GmbH De-Trey-Straße 1 78467 Konstanz Germany	 Maillefer Instruments Holding Sàrl Chemin du Verger 3 1338 Ballaigues Switzerland
	UK REP DENTSPLY IH Limited Brunel Way Stonehouse Gloucestershire GL10 3GB	 dentsplysirona.com/ifu
EU Importer		
Dentsply Sirona Europe GmbH Sirona Straße 1 5071 Wals bei Salzburg Austria		

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