## Grandio®SO Flow

### Grandio®SO Heavy Flow

In accordance with DIN EN ISO 4049

#### Instructions for use

##### Product description:

GrandioSO Flow / Heavy Flow is a flowable composite. In contrast to conventional flows, it is characterised by its increased filler content and stability. The radiopaque restorative contains 81 % (GrandioSO Flow), 82 % (GrandioSO Flow WO) and 83 % w/w (GrandioSO Heavy Flow) inorganic fillers in a methacrylate matrix and cures under blue light. In GrandioSO Flow WO, an extremely radiopaque version is at your disposal. GrandioSO Flow / Heavy Flow is used with a dentine-enamel bond.

##### Shades:

A1, A2, A3, A3.5, A4, B1, B2, C2, OA2, BL, WO, Incisal, GA5

##### Indications:

- Fillings of class I - V cavities
- Minimally invasive filling therapy
- Extended fissure sealing
- Blocking out undercuts
- Lining or coating cavities
- Repairing fillings and veneers
- Luting translucent prosthetic pieces (e.g., full ceramic crowns, etc.)
- Interlocking and splinting teeth (e. g., with glass fibre strands)
- Use as a base material in combination with glass fibres or similar products for fabricating semi-permanent crowns and bridges

##### Preparation:

Clean the teeth to be treated with a fluoride-free cleaning paste.

##### Shade selection:

To select a shade, compare the enclosed shade guide with the cleaned tooth before anaesthesia. Moisten the shade sample and tooth and select the shade, preferably in daylight.

##### Creating a dry working field:

Ensure that the working field is sufficiently dry. Use of a rubber dam is recommended.

##### Cavity preparation:

In principle, cavity preparation should be carried out according to the rules of adhesive filling therapy and should be minimally invasive to conserve healthy tooth substance. Bevel the enamel margins on anterior. Pay attention to special preparation shapes for extended fissure sealing and deciduous teeth! Clean and dry the cavity afterwards. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here.

##### Lining:

A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

##### Matrix placement:

Take appropriate measures for areas where the applied GrandioSO Flow / Heavy Flow can flow out. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and insertion of the matrix.

##### Handling:

Allow the material to reach room temperature before use.

##### Bonding material:

GrandioSO Flow / Heavy Flow is used in the adhesive technique with a dentine-enamel bond. All light-curing bond materials can be used. Follow the respective instructions for use with regard to preparation (etching technique) and application.

##### Application of GrandioSO Flow / Heavy Flow:

**Syringe:** Lock the enclosed application tip on the syringe with a clockwise turn and directly apply GrandioSO Flow / Heavy Flow. For fillings with a thickness of over 2 mm, apply and cure in layers. Do not pull the plunger back on the syringe during or after use to guarantee the function of the non-dripping syringe.

**Caps:** Remove protective cap. Insert Caps into the opening of the dispenser (follow the instructions for use). Bend the metal cannula according to the indication into the desired direction. Apply GrandioSO Flow / Heavy Flow directly into the cavity by slightly and evenly pressing the levers of the dispenser together. Use each GrandioSO Flow / Heavy Caps for one patient only!

##### Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. The curing time is as follows when using an LED-/halogen light with a minimum output of 500 mW/cm<sup>2</sup>:

20 s for shades A1, A2, A3, A3.5, B1, B2, C2, BL and Incisal

40 s for shades A4, OA2, WO and GA5

Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced. Incomplete curing may lead to discoloration and discomfort.

##### Finishing:

The restoration can be finished and polished immediately after removing the matrices (e.g., fine or extra-fine diamond grinder, polisher). The tooth should be fluoridated as a final step.

##### Information / precautionary measures:

- No known side effects. Sensitization in hypersensitive patients, however, cannot be excluded.
- GrandioSO Flow / Heavy Flow contains methacrylates and BHT. GrandioSO Flow / Heavy Flow should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- The enclosed application tips are intended for a single use and must be disposed of after only one use. Return the cap to the syringe after use.
- Phenolic substances, especially preparations containing eugenol or thymol, interfere with the curing of filling composites. Avoid the use of zinc oxide eugenol cements or other eugenol-containing materials in combination with filling composites.

##### Storage:

Store at temperatures between 4 °C - 23 °C. To avoid exposure to light and potentially resulting polymerisation, syringes should be closed immediately after dispensing. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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