- Information for User (Product Leaflet) -

PRUDENTI 
Innovative dental products

Pocket-X[®] gel

Information for User (Product Leaflet)

For use by the Dental Health Care Professional

Read this leaflet in its entirety before using this medical device.

Pocket-X Gel is a medical device in the form of a periodontal in-situ gelling sustained release subgingival delivery system. It is an adjunctive treatment after root surface debridement (RSD) or scaling, to improve the healing of the gingiva and to inhibit bacterial re-colonization in the periodontal pocket and as part of the routine periodontal treatment program performed by a dentist. Thanks to the mucoadhesive properties of the in-situ self-forming gel, it is formed and bulked in the application site and adheres to the mucosal tissues.

Mode of Action

Pocket-X Gel facilitates its effect via a physical mode of action. It is delivered into the periodontal pocket after scaling or root planing and cleaning of the pocket from plaque. Inside the pocket, the gel acts as a filler for volume restoration, thus physically preventing new bacteria from entering and re-infecting the pocket.

The gel is easily applied into the deepest portion of the periodontal pocket by means of a thin plastic tip or a thin flexible brush with a rounded tip that will not traumatize the tissues and ensure better application of the gel at the points which are deepest and most difficult to reach. Using Pocket-X Gel is easy due to simple physics; liquid fills pockets better than gel does, yet Pocket-X Gel is a liquid that transforms into a gel and stays at the affected site.

Composition:

Each syringe with 1 ml contains: Water, Poloxamer 407, Phenoxyethanol, Hyaluronic Acid, Octenidine Hydrochloride.

Intended use and indications:

Pocket-X Gel is intended to improve the healing of the gingiva and to inhibit bacterial re-colonization in the periodontal pocket. Pocket-X Gel is indicated to be used as an adjunctive treatment after root surface debridement (RSD) for moderate to severe chronic periodontal disease in adults with periodontal pockets. It is recommended also as part of the routine periodontal treatment program performed by a dentist.

Contraindications:

- Patients with known hypersensitivity to Hyaluronic Acid or to any of the ingredients listed in the composition section.
- Pocket-X Gel is not indicated for use in children and adolescents, since its safety and efficacy in this age group have not been established.

Special warnings and precautions for use:

Do not inject! Pocket-X Gel is intended for subgingival application only. Local hypersensitivity reactions such as gingival swelling might occur. Patients should be instructed to seek immediate medical attention if after exposure to the product they develop allergic symptoms such as: skin rash, itching, generalized swelling, breathing difficulties, light-headedness, rapid heart rate, upset stomach or diarrhoea, etc.

Interactions with other products and other forms of interactions:

There are no known interactions when administering Pocket-X Gel with other oral products or oral treatments. From review of the available literature, it is unlikely that problems will arise from using any other oral products together with Pocket-X Gel. However, it is important that the treating physician is informed and shall carefully consider the use of Pocket-X Gel in combination with other products.

As a precaution: Pocket-X Gel should not be used together with other products intended to act on the gums, teeth and in the oral cavity without prior consideration.

Fertility, pregnancy and lactation:

Pregnancy: The safety of Pocket-X Gel in pregnant women has not been established. No controlled clinical trials with Pocket-X Gel have been conducted

in pregnant women or in animals.

As a precaution: Pocket-X Gel should not be used in women who are pregnant, trying to become pregnant or suspect they are pregnant; and the administration of Pocket-X Gel should be suspended for the duration of pregnancy or until it has been determined that the woman is not pregnant. If administered, the dentist should weigh the expected benefits of the product to the mother against the possible risks to the foetus.

Lactation: There are no adequate experimental clinical and animal data on the use of Pocket-X Gel during lactation. As a precaution, women administered with Pocket-X Gel must not breast-feed their infants for a period of 3 weeks after the last administration of Pocket-X Gel.

Fertility: There are no human and animal data regarding the effect of Pocket-X Gel on fertility. As a precaution, Pocket-X Gel should not be used in patients trying to conceive. Fertile men and women should use a contraceptive for 3 weeks if Pocket-X Gel is to be administered.

Effects on ability to drive and use machines:

Pocket-X Gel is not likely to affect the ability to drive or use machines.

Posology and method of administration:

For use by dentists.

Dosage – adults including the elderly: following mechanical debridement (RSD or Scaling), the liquid is applied to each of the periodontal pockets to be treated and transforms within a few seconds into a gel filling the pocket. If needed, additional treatment can be performed, following additional mechanical debridement, within no less than 2 weeks after the previous treatment. Prior to administration of the gel, attach to the syringe the preferred dedicated tip that is supplied with the product. There is no clinical difference with which tip you choose to use. Both tips are suitable for use in the same periodontal pockets. The difference is only on the way each tip delivers the gel inside the periodontal pocket. You can switch between the tips to achieve the optimal ease of administration:

Option A: Red coloured tip with a small brush.

This is a short length tip with a small brush on its end, which helps to spread the gel in the periodontal pocket.

Option B: Blue coloured tip.

This is a long length tip, which enables easier placing of the gel in the inner (deepest) areas of the periodontal pocket.

Each syringe and tip are intended for use on a single patient in one or numerous periodontal pockets. Dispose of the syringe and tips after use according to local regulations. Do not use any gel remaining in a syringe on a different patient.

Pocket-X Gel contains enough volume for several pockets so it can be applied simultaneously to all the periodontal pockets that need to be treated. Do not exceed one full syringe per patient per treatment. The patient should be instructed to continue normal oral hygiene procedures. No restrictions on dietary habits are needed.


Dosage – children and adolescents: the gel is not indicated in children and adolescents, since safety and efficacy in these age groups have not been established.

Method of administration:

- 1) Remove one syringe from its package and unscrew its cap. Then screw to the syringe the preferred plastic tip (red coloured tip or blue coloured tip).
- 2) Hold the syringe and place the end of the tip on the top part of the periodontal pocket. Insert the tip into the periodontal pocket to its maximum depth.
- 3) Press slowly the syringe plunger so the liquid is released at the base of the pocket until the gel becomes visible at the gingival margin.
- 4) Remove the syringe from the mouth.

The gel's consistency allows the placement into the pocket with little discomfort to the patient. Pocket-X Gel is provided as a liquid pre-filled formula in a syringe

and transforms into a gel at body temperature. If it has formed into a gel prior to its application, e.g., due to warm room temperature, then the syringe should be placed in a refrigerator for a few minutes until it becomes liquid again. Notice, that the product's storage conditions state "Do not refrigerate!" yet you can safely place the product for a few minutes in the refrigerator in order to enable the gel to become liquid again. This does not affect the product's safety, efficacy or quality.

Instructions for the patients:

Do not eat for 2 hours after treatment. Do not rinse or gargle while drinking for 2 hours after treatment. Do not drink hot beverages for 4 hours after treatment. Brush the teeth very gently in the two days after treatment. Do not use accessories for cleaning between the teeth (e.g., dental floss, etc.) in the first 2 days after treatment.

Undesirable effects:

Some of the patients experience adverse reactions, usually transient, during the first few hours after application of the gel in the periodontal pocket. These may also be due to the mechanical application in the periodontal pocket or as a result of the preceding scaling procedures. The most common undesirable effects are discomfort in the treated area, upper gastrointestinal system disorders in the mouth, and dental, gingival, or oral soft tissue reactions.

Oral disorders: Toothache, gingival swelling, gingival pain, gingival bleeding, gingival hyperplasia, gingival recession, gingival pruritus, mouth ulceration, sensitivity of teeth.

Other undesirable effects may include skin disorders such as: Hypersensitivity to the product causing allergic skin reactions such as dermatitis, pruritus, erythema, eczema, rash, urticaria, skin irritation and blisters.

Other adverse reactions:

In rare cases patients may experience systemic hypersensitivity, anaphylactoid reaction, cellulitis and abscess related to the application site, loss of taste and gingival discoloration.

Reporting of suspected adverse reactions:

Reporting suspected adverse reactions is important. It allows continued monitoring of the benefit / risk balance of the products. Healthcare professionals are asked to report any suspected adverse reactions via the manufacturer details in this information leaflet, via the local distributor or via the European and/or Swiss Authorized Representative. Adverse events and suspected adverse events should also be reported to the local competent authority.

Overdose:

No case of overdose has been reported.

Shelf life:

The expiry date of the product is indicated on the syringe and on the box. Do not use the product after the expiry date stated on the package.

Shelf life after first opening the immediate packaging:

Use immediately after opening and discard after use.

Special precautions for storage:

Store at room temperature (15°C - 25°C). Keep in its original package to protect from light. Avoid exposure to direct heat and sunlight. Do not freeze! Do not refrigerate!

Keep out of reach of children.

Special precautions for disposal:

Dispose of the used syringes and tips in accordance with local regulations concerning healthcare waste.

Manufactured exclusively for: Prudentix Ltd., Israel.

Tree of Life Pharma Ltd., 1 Bat Sheva Street, P.O.B. 799, Lod 7110604, Israel.

MedNet EC-REP GmbH, Borkstrasse 10, 48163 Münster, Germany.

Tao of Excellence GmbH, Vorstadt 26, 8200 Schaffhausen, Switzerland.

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Geistlich Pharma AG, Bahnhofstrasse 40, 6110 Wolhusen, Switzerland.

| | | | |
|---------------|---|--|--|
| Symbols used: | | | Keep out of reach of children |
| | Consult instructions for use | | Batch number |
| | Caution | | Use-by date (expiry date) |
| | Do not re-use | | Medical Device |
| | Not to be swallowed! | | CE mark medical device certified by notified body number |
| | Not for injection! | | Manufacturer |
| | Maximum and minimum temperatures the device can be exposed to | | Authorised representative in the European Community / European Union |
| | Keep away from sunlight | | Authorised representative in Switzerland |
| | Do not freeze! | | Importer |
| | Do not refrigerate! | | Distributor |

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Tillverkad exklusivt för: Prudentix Ltd., Israel.
EC REP MedNet EC-REP GmbH, Borkstrasse 10, 48163 Münster, Tyskland.
CH REP Tao of Excellence GmbH, Vorstadt 26, 8200 Schaffhausen, Schweiz.
Geistlich Pharma AG, Bahnhofstrasse 40, 6110 Wolhusen, Schweiz.