

## **TePe Interspace**<sup>™</sup>

Handle Polypropylene (PP).

Filaments Polyamide (PA).

Neck Can be angled after heating.

Autoclave Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min,

sterilization pressure: 3 bar / 300 kPa.

Pack 1 handle + 12 tips/blister pack or 1 handle + 4 tips/blister pack. 12 blister packs/carton. The plastic front is made

of Polyethylene terphtalate (PET). The backing card is made of cardboard. The carton is made of chlorine-free

corrugated cardboard.

Recycling The brush in the usual household waste. The blister pack should be separated and recycled:

the plastic front in container for plastic waste and the backing card in container for recycling of cardboard.

The carton should be placed in container for corrugated cardboard.

**Product** - Angled brush with a pointed tip.

features - Interchangeable tips.

- Two textures: soft (yellow tip), medium (blue tip)

**Examples of** - Implants.

**areas of use** - Orthodontics appliances.

- Difficult-to-reach areas, e.g. furcations and the back tooth's

distal surface

**User tip** For optimal access two different angles can be chosen.

