VDW® Endo Easy Efficient

RECIPROC® Minima

EN	Instructions For Use	2
FR	Instructions d'utilisation	10
IT	Istruzioni per l'uso	18
DE	Gebrauchsanweisung	26
ES	Instrucciones de uso	34
PT	Instruções de utilização	42
NL	Gebruiksinstructies	50
SV	Användarinstruktioner	58
DA	Brugsanvisning	66
NO	Brukerinstruksjoner	74
FI	Käyttöohjeet	82
PL	Instrukcja obsługi	90
LT	Naudojimo instrukcija	98
LV	Lietošanas norādījumi	106
ET	Kasutusjuhised	114
CS	Návod k použití	122
SK	Návod na použitie	130
HU	Használati utasítás	138
RO	Instrucțiuni de utilizare	146
SR	Uputstva za upotrebu	154
HR	Upute za uporabu	162
BG	Инструкции за употреба	170
BS	Uputstva za upotrebu	179
MK	Упатство за употреба	187
EL	Οδηγίες χρήσης	195
UK	Інструкції із застосування	203
TR	Kullanım Kılavuzu	211
KO	사용 지침	219
ID	Petunjuk Penggunaan	226
PT-BR	Instruções de uso	234

VDW®
Endo
Easy
Efficient

RECIPROC® Minima

Instructions For Use

dentsplysirona.com/ifu

FOR DENTAL USE ONLY
FOR PROFESSIONAL USE ONLY
STERILE – STERILIZED BY RADIATION
SINGLE USE

INSTRUCTIONS FOR USE

SINGLE USE STERILE ENDODONTIC INSTRUMENTS

1 PRODUCT DESCRIPTION

- Product description:
 RECIPROC® Minima files are engine-driven endodontic instrument made of nickel-titanium alloy used to mechanically shape and prepare the root canals during endodontic therapy. These devices are sterile, single use root canal instruments.
- Clinical benefits:
 The clinical benefit for the patient being treated with RECIPROC® Minima is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Indication for use:
 The products are indicated for treatment of endodontic disease.
- Intended purpose:
 Engine-driven instrument intended for root canal preparation (shaping and debridement of the root canal).
- Intended patient population:
 Patients who need root canal treatment.
- Intended user and environment of use:
 Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

RECIPROC® Minima instruments are delivered in sterile blisters.

The RECIPROC® Minima range comprises the following endodontic instruments:

Code **xx-length** **yyy-file size** **z-quantity per blister**	**Product**	**ISO XXX YYYv** **XXX-ISO size** **YYY-taper size** **v-variable taper**	**Available lengths**	**Packaging and content**
M ST RCPM zxx yyy	RECIPROC® Minima	R20: ISO 020.008v R25: ISO 025.008v R35: ISO 035.006v R45: ISO 045.005v	21/25/31 mm	Blister with 6 instruments (sterile) Blister with 4 instruments (sterile)
M ST RCPM 6xx AST	RECIPROC® Minima Assortment	R20: ISO 020.008v R25: ISO 025.008v R35: ISO 035.006v R45: ISO 045.005v	21/25/31 mm	Blister with 2 instruments R20, 2 instruments R25, 1 instrument R35, 1 instrument R45
M ST RCPM 4xx AST	RECIPROC® Minima Assortment	R20: ISO 020.008v R25: ISO 025.008v R35: ISO 035.006v R45: ISO 045.005v	21/25/31 mm	Blister with 1 instrument of each instrument (sterile)
M DD RCPM 625 yyy	RECIPROC® Minima Demo Pack	R20: ISO 020.008v R25: ISO 025.008v R35: ISO 035.006v R45: ISO 045.005v	25 mm	Blister with 6 instruments (sterile)
M DD RCPM 4xx AST	RECIPROC® Minima Demo Pack Assortment	R20: ISO 020.008v R25: ISO 025.008v R35: ISO 035.006v R45: ISO 045.005v	21/25/31 mm	Blister with 1 instrument of each instrument (sterile)

5 COMPOSITION

- RECIPROC® Minima instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with colored rings.
- The files are manufactured with a process that results in a file that has a colored appearance. Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect, and it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

6 COMPATIBLE DEVICES

- RECIPROC® Minima instruments shall be used in combination with an endodontic motor with reciprocating motion. Only use a dedicated Dentsply Sirona motor and contra-angle with RECIPROC® settings (manually in severe curvatures).
- Lubricants such as NaOCl, EDTA, ProLube, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator is the appropriate method of working length determination. The stopper should be correctly positioned, with the aid of pliers on the occlusal / incisal reference point for an accurate length measurement.

7 SAFETY NOTES

7.1 Warnings

- Strictly follow the Instructions For Use for Endodontic Products to minimize the following risks to the device, the patient and/or the user:
 - Breakage of instrument.
 - Cross-contamination.
 - Heat generation due to insufficient lubrication and irrigation.
 - Swallowing of working part of the instrument.
 - Over instrumentation.
- Sterility cannot be guaranteed if the packaging is open, damaged or wet.
- Similar to all machine-driven root canal instruments, RECIPROC® Minima instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Visually inspect the integrity of the packaging (sterile barrier when applicable) before use. Do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Pay particular attention to cutting edges in order to avoid injury and damage to the products.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Use Working Length tools (e.g. ENDO-M-BLOC®, or MINI ENDO-BLOC) to properly position the silicone stop at the targeted working length.
- The access cavity shall always be filled with NaOCl.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution at a concentration not greater than 5%.
- Clean and inspect the instruments frequently during instrumentation (i.e. after each wave), and discard any instrument with sign of defects showing untwisting of the blades, elongation, or wear such as uneven flutes or damaged cutting edges. For challenging canals, it is recommended to clean the instruments more often to avoid the accumulation of debris in the flutes.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.

- Irrigate abundantly and frequently the canal throughout the procedure (according to good dental practices). We recommend the use of DS irrigation needle and/or SmartLite Pro EndoActivator®.
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
 - Use a new file and discard it after the canal was treated (single canal use).
 - Use manual instead of rotary files.
 - Use small size, flexible and/or NiTi files (this will help avoid canal transportation).
 - Visually inspect the working part after each wave and discard any instrument with sign of defects showing untwisting of the blades, elongation, or wear such as uneven flutes or damaged cutting edges.
 - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

8 STEP BY STEP INSTRUCTIONS

RECIPROC® Minima instruments are used with slow in- and out-movements (pecks). The amplitude of each peck must not exceed 3 mm. Only very light pressure should be applied; never force the instruments down the canal.

For narrow and challenging canals, the procedure should start with instrument size R20.

In cases of wider canal anatomy, the selection of the initial instrument size is at the clinician’s discretion.

1) Review different horizontally angulated radiographs, or CBCT images (if required), to assess the canal anatomy, width, and curvature, and estimate length.
2) Prepare an adequate access cavity to allow easy identification of the canal orifices.
3) Place NaOCl or EDTA gel (if the canals are narrow / calcified) in the access cavity.
4) Select R20 and set the stopper at 2/3 of the estimated length, which is determined from the radiographs or CBCT images.

5) Introduce R20 into the access cavity. Press the ON/OFF button or foot pedal of the motor when the instrument is at the canal orifice. Move the instrument in a slow in-and-out pecking motion. The amplitude of each peck must not exceed 3 mm. Only very light pressure shall be applied - never force the instrument down the canal. After 3 pecks, remove, clean and inspect R20.
6) Irrigate the canal, and confirm patency with an ISO size 010 or 008 hand instrument (K-File or C-PILOT®) taken to 1-2 mm beyond the length reached by R20.
7) Repeat steps # 5 and 6 until R20 reaches 2/3 of the canal (indicated by the silicone stopper), usually after 2-3 series of pecks. Remove the instrument.
8) Repeat steps # 5 and 6 until R20 reaches length indicated by the integrated apex locator. Otherwise, use an ISO size 010 or 008 hand instrument to determine length with a tabletop apex locator. Then, repeat steps # 5 and 6 until R20 reaches working length (indicated by the silicone stopper). Re-confirm working length electronically with R20.
9) If R20 stops progressing, do not force the instrument; remove it from the canal. A glide path should be created (e.g. with R-Pilot® instrument).
 a) Use an ISO size 008 or 010 instrument, with the irrigant in the canal, to obtain patency to length as determined and confirmed with an apex locator.
 b) Create a glide path (e.g. with R-Pilot®) with a dedicated DS motor. Use very gentle and pressure-less pecking. The amplitude of each peck must not exceed 3 mm. After 3 pecks, remove the glide path instrument, inspect and clean it, and irrigate. Continue in the same manner to working length (as indicated by the stopper).
 c) Once a glide path is established, complete preparation to working length with R20 by repeating steps # 5 and 6. Re-confirm working length with R20.
10) If required, apical preparation with a larger ISO size RECIPROC® Minima can be done. Select the appropriate instrument size for the canal.
11) Always proceed with 3-D disinfection protocols. We recommend using the SmartLite Pro EndoActivator®.
12) Dry the canal with dedicated paper points before obturation with dedicated matched gutta-percha points and sealer (we recommend using AH Plus sealer).

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- RECIPROC® Minima instruments are for single use only and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization are not applicable.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

VDW®
Endo
Easy
Efficient

CE 2797

Authorised representatives

EU REP
DENTSPLY DeTrey GmbH
De-Trey-Straße 1
78467 Konstanz
Germany

UK REP
DENTSPLY IH Limited
Brunel Way
Stonehouse
Gloucestershire
GL10 3GB

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