


<b>PRODUCT</b>	Hawe Matrices Micro Thin																	
<b>REFERENCE</b>	<table border="1"> <thead> <tr> <th colspan="3">Hawe Matrices Micro Thin</th> </tr> </thead> <tbody> <tr> <td rowspan="6">Refill pack of 30</td> <td>Hawe Matrices Micro Thin</td> <td>405</td> </tr> <tr> <td>Hawe Matrices Micro Thin</td> <td>406</td> </tr> <tr> <td>Hawe Matrices Micro Thin</td> <td>407</td> </tr> <tr> <td>Hawe Matrices Micro Thin</td> <td>408</td> </tr> <tr> <td>Hawe Matrices Micro Thin</td> <td>409</td> </tr> <tr> <td>Hawe Matrices Micro Thin</td> <td>410</td> </tr> </tbody> </table>			Hawe Matrices Micro Thin			Refill pack of 30	Hawe Matrices Micro Thin	405	Hawe Matrices Micro Thin	406	Hawe Matrices Micro Thin	407	Hawe Matrices Micro Thin	408	Hawe Matrices Micro Thin	409	Hawe Matrices Micro Thin
Hawe Matrices Micro Thin																		
Refill pack of 30	Hawe Matrices Micro Thin	405																
	Hawe Matrices Micro Thin	406																
	Hawe Matrices Micro Thin	407																
	Hawe Matrices Micro Thin	408																
	Hawe Matrices Micro Thin	409																
	Hawe Matrices Micro Thin	410																
<b>COMPANY/ Manufacturer</b>	KerrHawe SA																	
<b>CLASSIFICATION</b>	CLASS I MEDICAL DEVICE																	
<b>PRODUCT DESCRIPTION</b>	Stainless steel matrices																	
<b>SPECIFICATIONS</b>	<p>405:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> <li>Overall dimensions: 50.2 mm x 10.4 mm</li> <li>Thickness: 0.03 mm</li> </ul> <p>406:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> <li>Overall dimensions: 50.7 mm x 11 mm</li> <li>Thickness: 0.03 mm</li> </ul> <p>407:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> <li>Overall dimensions: 50.5 mm x 10.5 mm</li> <li>Thickness: 0.03 mm</li> </ul> <p>408:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> <li>Overall dimensions: 64.5 mm x 10.5 mm</li> <li>Thickness: 0.03 mm</li> </ul> <p>409:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> <li>Overall dimensions: 64.4 mm x 14.5 mm</li> <li>Thickness: 0.03 mm</li> </ul> <p>410:</p> <ul style="list-style-type: none"> <li>Material: Stainless steel 1.4310</li> </ul>																	

	<ul style="list-style-type: none"><li>• Overall dimensions: 66.8 mm x 11.1 mm</li><li>• Thickness: 0.03 mm</li></ul>
<b>PACKAGING</b>	Refill pack: 5.8 cm width 8.6 cm length, 2.9 cm height
<b>SHELF LIFE</b>	1800 days
<b>INTENDED OF USE</b>	Use for both molar and premolar restorations
<b>SECURITY</b>	Only professional use
<b>TRANSPORT AND STORAGE CONDITIONS</b>	Ambient temperature