


AH Plus® Bioceramic Sealer

What sets this sealer apart?

Bioactivity of bioceramic sealers is important but you need more to ensure a stable seal. The apical seal is only completed after the sealer sets.

AH Plus Bioceramic Sealer sets apart from the rest because it sets quickly, has a high washout resistance and is more radiopaque compared to current market leaders.



1. data on file, compared to the most common bioceramic sealers in the market.



AH Plus[®] Bioceramic Sealer

Safe and gap-free healing

AH Plus Bioceramic Sealer is a tricalcium silicate-based sealer that stimulate formation of hydroxylapatite that promote tissue healing in the root canal. The sealer also has the feature of being both lipophilic and hydrophilic which create a gap-free bond with zero shrinkage.

Dimensionally stable & high washout resistance

Thanks to its low solubility and low film thickness AH Plus Bioceramic Sealer stably seals the root canal for a sufficiant period of time.





Better radiopacity

AH Plus Bioceramic Sealer

Competitor bioceramic sealer

3. data on file.

25% more radiopaque than competitor bioceramic sealer brands¹.

Film Thickness³

25um

75um




1. data on file

Retreatability

Easily removable after setting with a hand file or Niti file.

Product reference

Starter Kit (3gr syringe/20 cannulas)	B KT AHBS 000 000
Syringe Refill (3gr syringe)	B 00 AHBS 100 03G
Cannula Refill (20 pc)	B 00 AHBS 020 24G





THE DENTAL SOLUTIONS COMPANY™