

# **AH** Temp<sup>TM</sup>

Calcium Hydroxide-Based Root Canal Dressing

 Direct and Precise Application into the Root Canal

- Reduced Desiccation
- Proven Compatibility with AH Plus

For better dentistry

DENSPLY

AH\_Plus\_folder.indd 1 2013-05-08 11:51

# AH Temp<sup>™</sup> combines the effective disinfection of root canals with easy handling due to the convenient AH Temp<sup>™</sup> Cannulas.

# The aim of Endodontic Treatments: Effective Disinfection of Root Canals


AH Temp™ is a highly alkaline paste with a pH level of 12 based on calcium hydroxide which is known for its antibacterial effect.

# **Reduced Desiccation**



- ✓ Less blocking of syringes over time
- Less waste of desiccated syringes

## Stable Consistency during Chairside Use<sup>1</sup>




Once the syringe cap is removed, AH Temp<sup>TM</sup> keeps its pasty consistency for 6 hours – no time pressure during the treatment.

In the daily routine, AH Temp™ is usable for 3 months



# Excellent Radiopacity enables reliable X-ray Control<sup>1</sup>



Data on file

AH\_Plus\_folder.indd 2 2013-05-08 11:51

<sup>&</sup>lt;sup>2</sup> Not a registered trademark of DENTSPLY International, Inc.

# Direct and Precise Application into the Root Canal Angulation 45° Working length 25 mm "The application of AH Temp™ is brilliant – it is the best root canal dressing available at the moment. In addition, it is supplied with an excellent cannula!" Dr. Matthias J. Roggendorf, University of Marburg, Germany External diameter 0.41 mm AH Temp™ Cannula Calxyl<sup>2</sup> Cannula Proven Compatibility with AH Plus Consistent High Bond Strength of AH Plus on Dentin<sup>1</sup> AH Temp™ is compatible with AH Plus® - the ideal combination of root canal dressing and definite root canal sealer. Proven consistent performance of AH Plus® after the removal of AH Temp $^{TM}$ from the root canal. AH Plus® without AH Plus® with previous use of previous use of $AH \; \mathsf{Temp}^{\mathsf{TM}}$ AH Temp™

AH\_Plus\_folder.indd 3 2013-05-08 11:51

 $<sup>^1</sup>$  M. Roggendorf, 2012, proven for application of AH Temp $^{\rm IM}$  and AH Plus $^{\rm B}$  Sealer, data on file.  $^2$  Not a registered trademark of DENTSPLY International, Inc.

# **Your Handling Benefits**







The No. 1 Endo Sealer<sup>1</sup> in a syringe with convenient application tips





#### Extra fast

- Increased efficiency no manual mixing anymore
- Better working ergonomics direct, precise placement



#### Extra easy

- Pre-filled, ready to use syringe
- Simple application, no assistance needed



#### Extra safe

- Clean application
- Precise placement of paste in the root canal

# **ReOrder Information**

#### AH Temp™ Kit

4 x AH Temp™ Syringes, 0.75 ml each

20 x AH Temp<sup>TM</sup> Cannulas

1x Finger Grip

1x Illustrated Technique Guide

1x Instructions For Use

ReOrder 606.20.123



## AH Plus Jet® Starter Kit

1 x AH Plus Jet® Mixing Syringe, 15 g

20 x Mixing Tips

20 x Intra Oral Tips

1 x Mixing Pad

1 x Illustrated Technique Guide

1 x Organizer Tray



20 x AH Temp™ Cannulas

ReOrder 606.20.124



## ReOrder 606.20.115

### AH Plus Jet®

Mixing Tips/Intra Oral Tips: 40 x Mixing Tips with Intra Oral Tip



#### ReOrder 606.20.116

#### AH Plus Jet® Refill

2 x AH Plus Jet® Mixing Syringe, 15 g

ReOrder 606.20.118



# **Disposa Shield**

500 x disposable barrier sleeves

ReOrder A8801351



#### **DENTSPLY Nordic Baltic**

Box 2024, Flygfältsgatan 12 SE-128 21 Skarpnäck Fax +46 (0) 8 447 60 55 www.dentsply.se

Phone, SE +46 (0) 8 685 65 05 Phone, DK +45 20 46 56 80 Phone, NO +47 976 08 438 Phone, FI +358 (0) 40 741 90 84 Phone, LV +371 29 45 59 45

For better dentistry



1 Data on file.