

Product information

SENDOLINE AB

S-File

The S cross section with double edges give the S-File excellent cutting ability.

It is designed with a thin cylindrical core and a non-cutting tip with a 90° angle and a 90° cutting angle. This design increases the flexibility and safety.

You can use the S-File with two different techniques: either “quarter turn” or usual filing.

The S-File is also available in nickel titanium named [NiTi S-File](#).

ISO/Sizes: 06-140

Length: 21 mm, 25 mm, 28 mm, 31 mm

Package: Individual or assorted (in sizes 15-40 or 45-80) in blister. All packages come with 6 pcs/blister

Material: Operating part: Stainless steel. Handle: Plastic. Stopper: Silicone.

