

TempoCemNE

- Optimal adhesion and removability
- Thin film thickness
- Automatic mixing and direct application



Temporary luting of provisional and permanent crowns, bridges, inlays and onlays

TempoCemNE and Luxatemp: the dream pair for temporary restorations. This temporary luting cement is free from eugenol – and can therefore be used for all temporary insertions, whether the final luting is conventional or adhesive. Thanks to its very thin film thickness (under $20\mu m!$), TempoCemNE provides a precise fit.

TempoCem is also available with eugenol and as TempoCem-Soft specifically for the cementation of high friction multi-part bridges.



Seating a temporary with TempoCemNE.



Removal of excess material after setting.

Package

Accessories		
1	Automix-Dispenser Type 25 1:1	REF 110253
50	Automix-Tips Blue 1:1	REF 110266
50	Smartmix-Tips Short	REF 212040
100	Automix-Tips Blue 1:1	REF 110267
Package		
	TempoCemNE-Automix	
1	Cartridge @ 60 g	
35	Automix-Tips	REF 110272
	TempoCemNE-Smartmix	
2	Smartmix-syringes @ 11 g	
20	Smartmix-Tips	REF 212104
	Package TempoCemNE-Handmix	
1	Tube @ 85 g base	
1	Tube @ 25 g catalyst	REF 110137

